



HN 763 Preformed Tube, Ivory PVC, North Nasal, Cuff Maxillo-facial tube

Sizes & Technical Data

- ☞ Anatomically shaped smooth preformed design gently curves the tube towards the patient's forehead, where it can be conveniently secured and connected to the breathing system
- ☞ The design ensures specific requirements of patients undergoing maxillo-facial surgery
- ☞ The preformed bend lies strategically close to nostrils and can be cut shorter for easier suctioning
- ☞ Manufactured from soft, ivory PVC to protect delicate mucosal tissues and facilitate intubation into the nasopharynx
- ☞ Clear markings at the levels of the vocal cords and nostrils assist correct placement
- ☞ Single use, no risk of cross infection
- ☞ The size of the tracheal tube internal diameter is printed on the Pilot Balloon.
- ☞ A black Intubation Depth Marker is located 3.0 cm proximal to the cuff. This assists the accurate placement of the tracheal tube tip within the trachea
- ☞ 15mm connector ensures full compatibility with circuit connections

Product Description	Product Code	I.D.mm	O.D. mm	Cuff resting diameter mm	Pack size
North Nasal Cuff Tube	HN 763-060	6.0	8.8	27	10
North Nasal Cuff Tube	HN 763-065	6.5	9.4	27	10
North Nasal Cuff Tube	HN 763-070	7.0	10.2	28	10
North Nasal Cuff Tube	HN 763-075	7.5	10.8	29	10
North Nasal Cuff Tube	HN 763-080	8.0	11.6	30	10

HANSRAJ NAYYAR MEDICAL INDIA

Administrative Office:

L8, Laxmi Industrial Estate, New Link Road,
Andheri (w), Mumbai - 400 053
Tel. # 022-26304581/82/83, 26346257
Fax # +91-22-26346258

Corporate Office:

1416, Maker Chamber V,
Nariman Point, Mumbai - 400 021
Tel: 022-22885375/, 22840912

Dehi Office:

30/29, Basement, East Patel Nagar,
New Delhi - 110 008
Tel: 011-42481530/31, Fax: +91-11 42481532

Chennai Office:

10, 2nd Floor, Raheja Complex,
Anna Salai, Chennai - 600 002
Tel: 044-28582501/2, Fax: +91-44 28582500

E-mail: mumbai@hansrajnayyar.com • Visit Us At: www.hansrajnayyar.com