



HN 713 Preformed Tubes, Clear PVC, South, Oral, Cuff
HN 710 Preformed Tubes, Clear PVC, South, Oral, Plain

Sizes & Technical Data

- ☞ An ideal tube for "South facing" intubations in Anaesthesia
- ☞ Clear tube allows moisture condensation from the breath to be seen clearly
- ☞ The Intubation length i.e. from tip to incisor teeth conforms to cut-to-length norms
- ☞ Smooth murphy eye avoids vocal cord damage during intubation
- ☞ Size identification on the "South" limb and on the sensitive pilot balloon
- ☞ Inflation line joins outside the mouth to avoid damage to teeth
- ☞ The contour of its limb avoids pressure on the chin and lower lip
- ☞ Its low-pressure cuff is designed to minimise tracheal trauma.
- ☞ Single use, no risk of cross infection
- ☞ The size of the tracheal tube internal diameter is printed on the Pilot Balloon.
- ☞ A black Intubation Depth Marker is located 3.0 cm proximal to the cuff. This assists the accurate placement of the tracheal tube tip within the trachea
- ☞ 15mm connector ensures full compatibility with circuit connections

Product Description	Product Code	I.D.mm	O.D. mm	Cuff resting diameter mm	Pack Size
South Oral Cuff Tube	HN 713-050	5.0	6.8	20	10
South Oral Cuff Tube	HN 713-055	5.5	7.4	20	10
South Oral Cuff Tube	HN 713-060	6.0	8.2	24	10
South Oral Cuff Tube	HN 713-065	6.5	8.8	24	10
South Oral Cuff Tube	HN 713-070	7.0	9.6	26	10
South Oral Cuff Tube	HN 713-075	7.5	10.2	26	10
South Oral Cuff Tube	HN 713-080	8.0	11.0	28	10
South Oral Cuff Tube	HN 713-085	8.5	11.6	28	10
South Oral Cuff Tube	HN 713-090	9.0	12.2	30	10
South Oral PlainTube	HN 710-030	3.0	4.2		10
South Oral PlainTube	HN 710-035	3.5	4.8		10
South Oral PlainTube	HN 710-040	4.0	5.4		10
South Oral PlainTube	HN 710-045	4.5	6.2		10
South Oral PlainTube	HN 710-050	5.0	6.8		10
South Oral PlainTube	HN 710-055	5.5	7.4		10
South Oral PlainTube	HN 710-060	6.0	8.2		10
South Oral PlainTube	HN 710-065	6.5	8.8		10
South Oral PlainTube	HN 710-070	7.0	9.6		10
South Oral PlainTube	HN 710-075	7.5	10.2		10